

by PolyMem®

SPORTSWRAP® ST

Benefits of SportsWrapST™ Roll

- Effective for use during activity
- Slimmer, thinner profile
- Helps manage strains, sprains, bruises and contusions where the skin is intact
- Adjunct to and synergistic with most physical modalities and other therapies
- Long lasting, drug-free pain relief



Simple and convenient to use.
Helps relieve pain, swelling and inflammation.

What experience have athletic trainers had with ST Rolls?

ST rolls were used successfully by athletic trainers of major league teams during baseball spring training 2010.



SportsWrapST™ Roll Benefits You and Your Patient

The SportsWrapST Roll provides the same benefits as the SportsWrap Rolls and Wraps, but in a slimmer, thinner profile for ease of use. SportsWrapST can be used for short-term response during activity and after injury or during initial recovery. ST helps reduce swelling, inflammation, bruising, stiffness or pain due to acute or chronic injury when the skin is intact. It is constructed of a patented polyurethane membrane material that can be used under supportive wrapping and taping with the same pain relieving benefits of SportsWrap green membrane.

Rigorous animal research studies show that SportsWrap, SportsWrapST and PolyMem dressings inhibit the action of nociceptors, the populous raw nerve endings found in the epidermis, thereby interrupting the mechanism of inflammation and pain.^{1,2}

Study results revealed that “There is robust, reproducible and statistically significant decrease in both secondary mechanical and thermal hyperalgesia” when wounds were wrapped with polymeric membrane dressings.¹ Human case studies have revealed dramatic reductions in edema, ecchymosis (bruising) and pain in a wide array of athletic injuries.³



Here’s How It’s Applied:

Strains, Sprains and Contusions – Where the Skin is Intact

Unroll the green membrane wrap in snug contact with the skin. Wrap in a continuous motion overlapping by 1/3 - 1/2. Only one layer is required under tape and wraps. Cut or tear when complete. Secure with self-adhesive tape, elastic wrap or garment. Remove following activity and showering. Discard the SportsWrapST and apply new.

For longer wear, apply the thicker, longer acting SportsWrap in order to obtain the same benefits during the day or at rest. If either SportsWrap or SportsWrapST is used overnight, remove in the morning and check the injury.




PRODUCT DESCRIPTION

REF

DIMENSIONS

PACKAGING



2.5" X 12' SportsWrapST Roll

9120

2.5" x 12' (6.3cm x 3.6m)

10 Per Bag
1 Bag Per Case

1. Alvin J. Beitz, PhD, Allison Newman, Alan R. Kahn, M.D., Timothy Ruggles & Laura Eikmejer. A Polymeric Membrane Dressing With Antinociceptive Properties: Analysis With a Rodent Model of Stab Wound Secondary Hyperalgesia; The Journal of Pain, February, 2004; 5(1):38-47.
2. Alan R. Kahn, M.D. A Superficial Cutaneous Dressing Inhibits Pain, Inflammation and Swelling In Deep Tissues; World Pain Conference, July 15-21, 2000.
3. Results obtained from human case study experience.



Ferris Mfg. Corp. 16W300 83rd Street, Burr Ridge, IL 60527 USA
Toll Free USA: 800.POLYMEM (765.9636) | International: +1 630.887.9797 | www.PolyMem.com

Unless otherwise indicated, all trademarks are owned by or licensed to Ferris.
© 2011 Ferris Mfg. Corp.

MKL-564, REV-0, 0211

MADE IN USA

