

Can a Drug-Free Dressing Decrease Inflammation and Wound Pain? What Does the Evidence Say?

Roger C. Sessions, DO, FACEP | Chairman and CEO, Ferris Mfg. Corp., 16W300 83rd Street, Burr Ridge, IL 60527 USA

PURPOSE

Persistent wound pain not only decreases the quality of life for the patient, but it also directly inhibits healing by producing a stress response, which promotes a catabolic state.^{1,2} In addition, wound pain predisposes the patient for wound infection through depression of the immune system.² Patients with multiple comorbidities and the elderly are often poor candidates for systemic pharmaceutical pain relievers. Even when they are prescribed, systemic products do not always address wound pain well.

Topical pharmaceutical pain relievers such as ibuprofen and lidocaine have local side-effects that can increase infection risks and/or directly impair wound healing.³⁻⁹ Studies indicate that inflammation often plays a critical role in the etiology of many wound types, but anti-inflammatory drugs such as prednisone inhibit wound healing. Can any intervention truly promote wound healing while decreasing inflammation and persistent wound pain?

RATIONALE

Pain and trauma associated with dressing changes (procedural pain) can be minimized by using a non-adherent dressing,¹⁰ such as polymeric membrane dressings.^{11,12} Components of these sophisticated dressings also work synergistically to provide continuous cleansing of the wound bed, usually completely eliminating the need for manual wound bed cleansing or even rinsing at dressing changes.¹²⁻¹⁶ Wound bed cleansing is an often overlooked source of procedural pain.¹⁷

But this study focuses on persistent wound pain, including examining the scientific physiological explanation for claims that polymeric membrane dressings inhibit the nociceptor response, even through intact skin.^{12-16,18} The resulting decreased pain and inflammation should dramatically improve wound healing.^{1,2} In addition, any drug-free persistent wound pain relief would avoid negative side-effects inherent in the use of currently available systemic and topical pharmaceutical agents.

METHODOLOGY

A thorough search of the peer-reviewed published data discussing this dressing's possible effect on inflammation and wound pain was performed to determine the robustness of the evidence for its ability to decrease inflammation and persistent wound pain. We identified four scientific studies: a 72-patient controlled study of pain and comfort level on burn and split-thickness skin graft patients, a 24-patient controlled study of pain and inflammation following arthroscopy and two rigorous animal studies performed by independent laboratories. We also examined the four published large facility-based product evaluations. Finally, we did a meta-analysis of every published peer-reviewed patient case study or series (187 individual patients) to identify how polymeric membrane dressings may influence persistent wound pain in clinical practice.

OBJECTIVES

1. Recognize that side effects of many pharmaceuticals currently available for wound pain can hinder healing.³⁻⁹
2. Review the evidence found in patient case studies and in rigorous scientific studies, for polymeric membrane dressings' ability to influence the nociceptor response, which can often decrease wound pain.
3. Consider the advantages of using a drug-free dressing over using systemic or topical medications for the relief of wound pain.

RESULTS

Laboratory research unequivocally concludes that polymeric membrane dressings achieve significant pain relief by inhibiting the nociceptor response both locally and centrally.¹⁹⁻²¹ It appears that the dressing is also influencing other receptors, resulting in additional analgesia beyond that anticipated by its antinociceptive properties (see data, below).²⁰⁻²¹

The facility-based evaluations found decreased pain, spasticity and bruising and increased patient mobility when polymeric membrane dressings were used. The decrease in persistent wound pain provided by polymeric membrane dressings was one of the critical factors influencing the choice to use polymeric membrane dressings (see map) in each of the four large evaluations.^{13,16, 23, 24}


Persistent wound pain was mentioned as a problem for 88 of the 187 case study/series patients. The independent clinician authors attributed a reduction in persistent wound pain to polymeric membrane dressings in 83 of these patients (94%). In 67 cases (76%), the pain reduction was described as dramatic, very significant, etc. (see chart, right).

CONCLUSION


Polymeric membrane dressings definitely inhibit the nerve conduction that normally leads to pain and the inflammatory response, not only on wounds, but also when applied to intact skin. This reduction in the nociceptor response occurs without interfering with the robust localized inflammatory response required for healing the injury. These insights into the actions of this dressing formulation have exciting implications for persistent wound pain relief and wound healing in both chronic and acute wounds.

CASE STUDY META-ANALYSIS

Pain Reduction: 83 Of those: Dramatic Pain Reduction: 67
No Data: 3 Same Pain: 1 Worse: 1




FACILITY-BASED EVALUATIONS



BURN AND GRAFT PATIENT STUDY¹¹

A Korean group investigated the effects of polymeric membrane dressings on wound healing and pain in 1999. Their work was published in a peer-reviewed journal, but since it was in Korean, it went largely unnoticed. Later, the authors directed a translation into English, from which this information was obtained. Using anesthetized rabbits, the authors determined that epithelialization was significantly increased with polymeric membrane dressings (p<0.05). They went on to test polymeric membrane dressings in a 72-patient controlled study of pain, healing time and comfort level on patients with burns and on patients with split-thickness skin graft (STSG) donor sites. Polymeric membrane dressings were significantly superior to standard care for each of the three parameters (p<0.01). The authors reported that the dressings were also easy to use and economical.




BLUNT TRAUMA ANIMAL MODEL¹⁹

A woman fell from a ladder, leaving a significantly painful area with a very slight abrasion. She placed a polymeric membrane dressing on the area, hoping that it would alleviate the pain. The woman experienced a complete absence of tenderness, swelling and bruising in the area covered by the dressing, but the injured area that was left exposed exhibited the expected response – note the line of ecchymosis (photo, right).

The surprising dressing-shaped area without ecchymosis exhibited by the woman who fell in the office, coupled with amazing results when polymeric membrane dressings were used post-arthroscopic knee surgery on a 65-year-old man, led Dr. Kahn of the University of Minnesota to do studies using the dressings over intact skin. The results were presented at the World Pain Conference in 2000.

- Uniform blunt trauma was applied to both legs on 14 anesthetized animals
- A polymeric membrane dressing plus a compression wrap was applied to one leg
- Only the compression wrap was applied to the other leg
- Two independent blinded observers evaluated swelling on a 0 – 4 scale at 24 and 48 hours


The author concluded: "Preliminary experimental evidence suggests that this cutaneous dressing inhibits the activity of the nociceptive neurons in the epithelium which, in turn, blocks the response of the spinal dorsal root mechanism that is responsible for generating swelling, inflammation and pain."¹⁹



ARTHROSCOPIC KNEE SURGERY²²

A 24-patient randomized controlled study of pain and inflammation using polymeric membrane dressings following arthroscopy was performed. The results were published in the Medline indexed journal, *Orthopedics*, in 2003.

Pain and skin temperature (an indicator of inflammation) were compared in 24 patients randomized to the use of polymeric membrane wraps or standard dressings following minor arthroscopic knee surgery over a ten day period. The treating surgeon, who performed the post-operative evaluation on day 10, was blinded. Despite equal use of pain medication, pain scores were significantly lower for the polymeric membrane group on each of the ten days, when compared to the controls. The polymeric membrane group demonstrated lower overall pain scores (p=0.03) and also lower temperatures (p=0.02) than the control group.



INCISIONAL STUDY ON A RODENT MODEL^{20,21}





Preliminary findings (graphs below) showed a statistically significant decrease in spinal cord Fos labeling when animals with incisions wrapped with polymeric membrane dressings were compared to those wrapped with a placebo foam¹ (p=0.005), gauze (p<0.0001) and no dressing at all (p=0.0042). The authors suggested that, while clearly the incised animals were exhibiting decreased nociceptor activity, polymeric membrane dressings must be activating non-nociceptive peripheral nerve fibers in the animals without incisions.

At the conclusion of the main study, the pain specialists hypothesized that polymeric membrane dressings use the same mechanisms as acupuncture and nitrous oxide administration. But, polymeric membrane dressings are non-invasive and do not have systemic effects.

¹the same foam substrate found in polymeric membrane dressings, but without the added components

Surprisingly, when polymeric membrane dressings were applied to animals without incisions in the preliminary study, they showed significantly more Fos labeling than animals wrapped in the placebo foam¹ (p<0.0001) or no dressing at all (p=0.0042). The authors suggested that, while clearly the incised animals were exhibiting decreased nociceptor activity, polymeric membrane dressings must be activating non-nociceptive peripheral nerve fibers in the animals without incisions.

At the conclusion of the main study, the pain specialists hypothesized that polymeric membrane dressings use the same mechanisms as acupuncture and nitrous oxide administration. But, polymeric membrane dressings are non-invasive and do not have systemic effects.



REFERENCES

1. Clay CS, Chen, WYJ. Wound pain: the need for a more understanding approach. *Journal of Wound Care*. 2005;14:4,181-184.
2. Middleton C. Understanding the physiological effects of unrelieved pain. *Nursing Times*. September 16, 2003; 99(37):28-31.
3. Divedi S, Tiwari SM, Shama A. Effect of ibuprofen and diclofenac sodium on experimental wound healing. *Indian Journal of Experimental Biology*, November 1997; 35:1243-1245
4. Proper SA, Fenske NA, Burnett SM, Luria LW. Compromised wound repair caused by perioperative use of ibuprofen. *Journal of the American Academy of Dermatology*. 1988; 18:1173-1179.
5. Brooks G, Yu X-M, Wang Y, Crabbe MJC, Shattock MJ, Harper JV. Non-steroidal anti-inflammatory drugs (NSAIDs) inhibit vascular smooth cell proliferation via differential effects on the cell cycle. *Journal of Pharmacy and Pharmacology*. 2003; 55:519-526.
6. Salcido R, et al. Do anti-inflammatories have a role in wound healing? *Adv Skin Wound Care*. 2005 Mar; 18(2):65-61. Session 301 Clinical Symposium on Advances in Skin and Wound Care, Oct. 23 – 26, 2005.
7. Jones MK, Wang H, Peskar BM, Levin E, Sarfeh IJ, Tamowski AS. Inhibition of angiogenesis by nonsteroidal anti-inflammatory drugs: insight into mechanisms and implications for cancer growth and ulcer healing. *Nat Med*. 1999 Dec; 5(12):1418-23
8. Tamowski AS, Jones MK. Inhibition of angiogenesis by NSAIDs: molecular mechanisms and clinical implications. *J Mol Med*. 2003; 81(10): p. 627-36.
9. Powell DM, Rodeheaver GT, Foresman PA, Hanks CL, Bellian KT, Zimmer CA, Becker DG, Edlich RF. Damage to tissue defenses by EMLA® cream. *Journal of Emergency Medicine*. 1999; 9:205-209.
10. European Wound Management Association (EWMA). Position Document: Pain at wound dressing changes. London:MEP Ltd, 2002.
11. Kim Y, Lee S, Hong S, Lee H, Kim E. The effects of polymem on the wound healing. *J Korean Soc Plast Reconstr Surg* 1999; 109:1165-1172.
12. Sessions RC. Severe road rash abrasions resulting in edema, inflammation, and pain in track race cycling athlete. *National Athletic Trainers Association (NATA) 58th Annual Meeting & Clinical Symposia*, June 26 – 30, 2007. Anaheim, CA, USA.
13. Tamir J. Polymeric foam dressing for skin grafts donor sites: 3 years experience on 800 cases. *23rd Clinical Symposium on Advances in Skin and Wound Care*. Poster #49, Oct 27 – 30, 2008. Las Vegas, NV, USA.
14. Sessions RC. Full-thickness chin wound healed in 14 days using only polymeric membrane dressings. *WOCN Society 39th Annual Conference*. Poster #1277, June 10 – 13, 2007. Salt Lake City, UT, USA.
15. Hubbard M. Pain relief and healing using polymeric membrane dressings under compression for venous hypertension ulcers. *20th Annual Symposium on Advanced Wound Care (SAWC)*. Poster #160, Apr 28 – May 1, 2007. Tampa, FL, USA.
16. Stenius M. Fast healing of pressure ulcers in spinal cord injured (SCI) people through the use of PolyMem® dressings. *10th anniversary EPUAP open meeting*. Poster #21, August 30 – September 1, 2007. Oxford, England.
17. Fleck, CA. Managing wound pain: today and in the future. *Advances in Skin and Wound Care* 2007;20:3,138-145.
18. Benskin L. Dramatic pain relief through the use of polymeric membrane dressings (with and without silver) on a deep axillary wound. *WOCN Society 38th Annual Conference*. Poster #167/Abstract #1686, June 24 – 28, 2006. Minneapolis, MN, USA.
19. Kahn AR. A Superficial Cutaneous Dressing Inhibits Pain, Inflammation and Swelling In Deep Tissues. *World Pain Conference*, July 15-21, 2000. Pain Medicine 2000 June;1(2):187.
20. Hayden JK, Cole BJ. The effectiveness of a pain wrap compared to a standard dressing on the reduction of postoperative morbidity following routine knee arthroscopy: a prospective randomized single-blind study. *Orthopedics*. 2003 Jan; 26(1):59-63.
21. Beltz A, Kahn A, Ferris PolyMem Plus™ dressing (REF 0548) Initial Study Summary. University of Minnesota. April 23, 2001. Unpublished.
22. Beltz AJ, Newman A, Kahn AR, Ruggles T, Elkmeier L. A polymeric membrane dressing with antinociceptive properties: analysis with a rodent model of stab wound secondary hyperalgesia. *J Pain*. 2004 Feb; 5(1):38-47.
23. Benskin L, Bolhuis J. Evidence-based skin tear protocol yields phenomenal results. *WOCN Society 40th Annual Conference*. Poster # /Abstract #2403, June 21-25, 2008. Orlando, FL, USA.
24. Wilson D. Skin tear healing improved through the use of polymeric membrane dressings. *21st Clinical Symposium on Advances in Skin & Wound Care*. Poster #341. 2006;373.

This research analysis was funded by Ferris Mfg. Corp.

¹PolyMem® dressings are made by Ferris Mfg. Corp., Burr Ridge, IL 60527 USA • 800/POLYMEM